

For the fourth NGS class, you will need to perform the following actions. As before, these links are also present on the course web site:

<http://physiology.med.cornell.edu/faculty/mason/lab/data/NGS/>

## 1. ChIP-Seq

Download the ChIPseeqer package:

<http://physiology.med.cornell.edu/faculty/elemento/lab/files/ChIPseeqer-1.0.zip>

Unzip the package and move into the directory

```
]unzip ChIPseeqer-1.0.zip  
]cd ChIPseeqer-1.0/
```

Move into the directory "libmd" and compile the library

```
]cd libmd/  
]perl Makefile.PL  
]make
```

Move back into the ChIPseeqer directory and compile chip-seeqer

```
]cd ..  
]make  
(OR, on Linux, type cd ../ ; make -f Makefile.linux)
```

Set environmental variable (make sure to use back-ticks! cmd substitiuon)  
export CHIPSEEQERDIR=`pwd`

Make a new directory for data analysis and move into that directory

```
]cd ..  
]mkdir CHIP  
]cd CHIP
```

Download raw data into the CHIP directory, taken from

[ftp://ftp.ncbi.nih.gov/pub/geo/DATA/supplementary/samples/GSM327nnn/GSM327662/GSM327662\\_hES.H3K4me3.aligned.txt.gz](ftp://ftp.ncbi.nih.gov/pub/geo/DATA/supplementary/samples/GSM327nnn/GSM327662/GSM327662_hES.H3K4me3.aligned.txt.gz)

Then unzip the chip data

```
]gunzip GSM327662_hES.H3K4me3.aligned.txt.gz
```

## 2. Prepare a genome assembler

Download Velvet

[http://www.ebi.ac.uk/~zerbino/velvet/velvet\\_0.7.62.tgz](http://www.ebi.ac.uk/~zerbino/velvet/velvet_0.7.62.tgz)

Gunzip and Untar the tarball

```
]gunzip velvet_0.7.62.tgz  
]tar -xvf velvet_0.7.62.tar
```

Compile the program

```
]make
```